

## 14in1 Hydra Beauty Skin System



### **Principle of Hydro Dermabrasion:**

Non-surgical skin refinish procedure, by using diamond dermabrasion heads to abrade or rub off the top skin layer, then vacuuming out the particles along with any dirt and dead skin. This procedure removes skin debris, imperfections, blemishes, wrinkles and unwanted pigmentation on the skin. Powerful diamond head microdermabrasion meets vacuum massage and intense aquation to provide gentle and effective exfoliation. Suitable for all skin types, including sensitive, rosacea and acne-prone skin.

# Instrument Functions



Water  
Dermabrasion



Skin Scrubber



Diamond  
Dermabrasion



Bubble pen



Spray Gun



Ultrasonic Face



Ultrasonic Eye



Bio



Tripolar RF



High Frequency



Vacuum Pen



Cold Hammer



PDT



Dr Pen



**Aqua peeling** The latest in skincare technology. Hydr dermabrasion utilizes the natural healing powers of water and oxygen to effortlessly exfoliate the skin without the use of hard crystals or abrasive textured wands, generating deeply hydrated and healthier looking skin..



### **Skin scrubber**

use of ultrasonic technology will change to 28000 times per second of electric shock waves to 28000 times per second of mechanical vibration wave the muscles of the body to me have osmosis.Can make skin more deep cells produce excited,so long as you use this instrument is used together with the corresponding potions or essence effect,can quickly eliminate shading wrinkles whelk corneous layer,able to pore dirt thoroughly clean skin,make skin shining white,recover the flexibility.



### **Diomand dermabrasion**

The clean and exfoliated skin is very receptive and absorbs maximally any treatments applied thereafter. Specific anti-oxidants,anti-aging, anti-blemish, anti-acne and many other specific treatments are then applied after the procedure for better penetration and results. This is not a simple facial, it is a professionally-administered, invigorating, advanced and medical-grade skin treatment that can change your skin forever



### **Bubble pen**

Magic bubbles take on an altogether new aspect of the oxygen, deep cleansing with rich oxygen foam, and deep cleansing of the skin, so that the skin can be whitening, refreshing and bright.




### **High Pressure Oxygen Spray**

The activity of oxygen and nutrients essential to promote through pulse-type deep into the skin ,effectively accelerated into cellular metabolism.It can make the lotion more delicate, so the skin will be easier to absorb, and effectively hyd rate,makes skin smoother.



### **Facial&eyes ultrasound**

Ultrasonic waves with deep penetration are produced by high-frequency oscillation of 1 million times per second. It is a kind of special instrument that emits a kind of waveform that can be transmitted by surrounding medium. It has more powerful energy than ordinary sound waves. It has high frequency, good directionality, strong penetrability and large expansion force.



### **Bio-polar**

Microcurrent is a safe and effective technology for clients who wish to attain/maintain a healthier and younger looking appearance.

Benefits :



### **Tripolar RF**

Non-invasive facelift, Reduce and eliminate fine lines and wrinkles  
Improve skin circulation, Improve hyperpigmentation and acne  
Increase product penetration



### **High frequency**

High frequency function is beauty salon acne-prone skin treatment good assistant of the high frequency

current heating can penetrate the skin fiber, promote the blood circulation, help skin absorb nutrition, eliminate waste, the heat can relax nerves, helping solvents penetrate the skin, reduce skin secretion, and sterilization, the role of promoting wound healing, effect is remarkable.

1. use a special ultraviolet radiation for local or punctiform light, have antiseptic effect, accelerate skin wound healing.
- 2 remove surface spots, acne and treat




### **Vacuum pen**

Vacuum negative pressure, in conjunction with the liquid to wash pores and peel skin surface, the black head and the cutin of skin away from the face to achieve a clean effect. Deep Cleaning, exfoliation, vacuum suction of all kinds of dirty in the skin.




### **Cold and Hot Hammer**

Smart temperature control system. Cold hammer is good for sooth and calm skin.



### **PDT light**

Photon dynamical nurse therapy is a new technology, using the photon power to do skin maintenance and treatment, transfer light energy to cell energy, then accelerate cell growth and blood circulation, stimulate fiber cells to produce collagen proteins, increase skin elasticity and take important part in rehabbing aging and acne skin, lightening flecks, skin upgrade and tighten and easing burnable skin..



### **Auto Micro Needle Pen**

Auto Micro Needle Pen uses multiple needles that vertically pierce the skin. This enhances rejuvenation results and is far safer for the client because there is much less epidermal damage. The automatic vibrating function of the Auto Micro Needle Pen increases the effectiveness of treatment by increasing absorption of products, while reducing pain and discomfort. Although non-ablative, Auto Micro Needle Pen is just as effective as treatments such as Fractional Laser Therapy, IPL, Laser Resurfacing, and Chemical Peels. The most amazing benefit is that Auto Micro Needle Pen's pricing is a fraction of the leading alternatives.

|                              |   |
|------------------------------|---|
| Product Name                 | 14 In 1 Hydra facial Machine  |
| Handles                      | 1. water dermabrasion with 6tips<br>2. Skin scrubber<br>3. Diamond dermabrasion with 9tips<br>4. Bubble pen<br>5. Spray gun<br>6. Ultrasonic face<br>7. Ultrasonic eyes<br>8. Bio<br>9. Tripolar RF<br>10. High frequency<br>11. Vacuum pen (3pcs)<br>12. Cold hammer<br>13. Photo light(PDT 2pcs)<br>14. Derma pen |
| Solution Bottle              | A, B, C,  |
| Bubble Bottle                | D   |
| Display                      | 15 inch Lcd touch screen  |
| The Power of the machine     | <350W   |
| Aqua peeling handle pressure | 90Kpa   |
| G.W                          | 45kg  |
| Packing size                 | 105*50*58cm Aluminum package  |

## Functions



# **HYDROGEN AND OXYGEN BUBBLE PEN SUCTION STRAIGHT TO DIG THE ROOTS**


DEEP CLEANS PORES, REMOVES CUTICLES AND  
ACNE BLACKHEADS, AND REMOVES OIL





**SPRAY GUN**  
**OXYGEN INJECTION AND HYDRATION**  
**REACH THE BOTTOM OF THE MUSCLE**


MOISTURIZES, DIMINISHES FINE LINES AND UNCLOGS PORES





**PHOTON LAMP**  
**TECHNOLOGICAL SKIN CARE**  
**IN JUST 3 STEPS**

TOPICAL USE, WHATEVER YOU WANT,  
MAKE EVERY INCH OF SKIN IMPECCABLE



Vacuum pen

It uses vacuum/suction to suck the blackheads from the pores, making our pores cleaner. And it's massaging during working, good for lymphatic drainage and active cells.



Promote blood circulation



Remove  
blackheads



Promote skin  
absorption



## Machine display



Right side



Left side



Back side



Production