



Alice 17 in 1 synthetic instrument

NV-WA16



01 Small bubble

Breaking through the technical level of traditional small bubble instruments, opening pores, softening keratin, cleansing the skin deeper, improving dullness, and blackheads leaving the skin clean, translucent, and hydrated.



02 V-Face

Using 20Hz-1000Hz frequency conversion EMS, the face is lifted and tightened, improving fine lines in the eyes, and restoring facial expression!




03 Five polar EMS

Using 20Hz-1000Hz frequency conversion EMS, the face is lifted and tightened, improving fine lines in the eyes, and restoring Eyes expression!



10 Radiant oxygen bubbles

Create an oxygen supply environment, activate cells to enhance absorption, and inject products into deep layers of the skin; Collaborate with fruit acid skincare products.



11 Beauty Spotted Pen

Remove and sweep spots, improve skin radiance

Lead-in

Recovery elasticity

Lifting

Beauty

Acne blackhead remove relieve auergy Deep
moistrization Deep clean tighten & nourish skin
moisturization Skin
rejuvenation increase elasticity



04 High pressure oxygen injection gun

High pressure oxygen injection penetrates deep into the skin layer, making it easier to absorb by narrowing pores.



05 Flat head Ultrasound

Using sound wave resonance frequency of 1-3 million vibrations to open pores and allow products to penetrate deep cells.



06 Bipolar warm wave

Used for facial thermal therapy to allow skin tissue to absorb a large amount of heat energy, combined with new collagen, to smooth and tighten the skin.



12 Ozone ion

Anti inflammatory, antibacterial, soothing and relieving itching, targeting allergic skin, improving skin condition



13 Diamond Dermabrasion

The principle of diamond particle friction and vacuum suction, Remove the aging keratin from the outer layer of the skin epidermis, To promote skin renewal, non-invasive Non irritating, safe, and convenient



14 Electric microcrystals

The front-end microcrystals or probes impact the surface of the skin back and forth, Generate channels for product import, Effectively promoting product absorption



Lifting



07 Point vibration shovel

The shovel handle utilizes high-frequency vibration mechanical waves to soften the skin's keratinocytes and deeply clean skin dirt.



08 Cold import

Quickly conduct skin qualitative analysis and achieve the effect of removing redness, anti allergy, and firming and shaping.



09 Quadrupole roller

Using 20Hz-1000Hz frequency conversion EMS, four roller massage, facial lifting and firming, improving fine lines in the eyes, and restoring facial radiance!




15 Point head Ultrasound

The front-end microcrystals or probes impact the surface of the skin back and forth, Generate channels for product import, Effectively promoting product absorption



16 LED colorful light mask

Red light promotes collagen regeneration and anti-aging
Blue light sterilization, repairing damaged skin
Green light improves redness and pigment redness
White light accelerates tissue metabolism and improves fine lines



17 High frequency wave

Anti inflammatory, bactericidal, soothing, and relieving itching
For allergic skin, Improve skin



For sensitive skin



For normal skin





For skin with large pores



Cleaning



Oily skin



Blackhead

Newfacial High end Water Drop Bubble

Specifications

Voltage: 100~240V 50Hz /60Hz

Power: ≤100W

Cold: 0-6°C

Ultrasound: 28~30KHz

HR Frequency: 0.2-0.5 Mhz

Ultrasound: 1.1MHz

Clean water mill: Positive pressure 4-5Kg
Negative pressure 90Kpa. Flow ≤30L

Display: 12.1 inches 1024x600

Small bubble product bottles: A, B, C, D
four types of bottles

TECHNOLOGY CREATE THE BEAUTY OF FUTURE

2024
AESTHETIC SERIES
Health & Beauty Creator



THANKS

