



## Triple Wave hair removal device



## THREE-BAND

HAIR REMOVAL + SKIN REJUVENATION

*Freezing point cooling*

*Rapid hair removal*

*Safe permanent painless*

*photorejuvenation*



## THREE-BAND DIODE HAIR REMOVAL

[755nm] [808nm] [1064nm]



808Diode laser hair removal

1000W 10Bar

High power




























## LASER PICOSECOND SOPRANO HAIR REMOVAL

### 2-in-1 Instruments

Eyebrow Washing  
Tattoo Cleaning  
Whole body available  
skin rejuvenation and firming

Exclusive Appearance

Your favor Our beauty mission







## Product parameters

Product name	808 Tri-Band Hair Removal + Picoseconds
Model number	DY T601
Package size	62*59*130cm
Weight of case	66kg
packing	Air cargo case
Hair Removal	755/808/1064nm
Laser Waveforms	532nm, 1064nm, 1320nm Non-invasive 532nm, non-invasive 1064nm
quantity	One machine

quantity	One machine
Screen size	10.4in
voltage	110V/220V
Cooling system	semiconductor cooling Air cooling + water cooling


## Function introduction

### Principle of operation

Through the epidermis to the hair follicle, according to the principle of selective light absorption, the energy of the laser is preferentially absorbed by the melanin in the hair, and then the hair loses its ability to regenerate, ensuring that the epidermis is not burned







### SHR STACK

Suitable for treatment of small areas such as the face (hairline, lip hair) bikini and fingers of the hands (feet)

### SHR

It is suitable for the treatment of large hairy areas. Characterized by: low energy density, high average power

### HR

For hairline, beard, bikini limbs and torso.



Suitable for many skin tones



I II III IV V VI

755nm for white skin (fine, golden hair removal )


808nm Optimum wavelength (Yellow, skin tone )

1064nm for black skin (black hair removal)



## Principle of Picosecond

The use of laser to generate a specific wavelength of intense pulsed light, will irradiate the tattoo area, it can penetrate into the deep layers of the skin, through the photothermal conversion, will be broken up by the pigment, so that with the body's metabolic activities to discharge, so as to achieve the effect of tattoo cleaning!





# 10000000+

Number of turns

Fast-breaking melanin



**532nm probe**

Skin spots, tattoo washing  
(red, brown, colored pigmentation)



**Non-invasive 532nm**

Non - invasive tattoo removal, for red  
coffee, color pigments, etc.



**Non-invasive 1064nm**

Non-invasive eyebrow wash  
for black blue, dark eyebrows



**1064nm**

Freckle removal, pigmentation, tattoo  
washing (black, blue, dark pigment)



**Carbon Laser/1320nm**


beautify white and tender skin, shrink pores,  
lighten spots and change skin quality







## Rejuvenation and Firming



Penetrates deep into the skin  
and improves skin texture





## CHIP UPGRADE

The use of high technology, the use of micro computer control perfect convenient and quick



## SAFETY GUARANTEE

One-button emergency stop button, instant stop machine, stable and safe



## LED TOUCH SCREEN

The picture is clearer and the whole is





more concise, the operation is more convenient, and the touch is more sensitive.

