





The image shows a white and black multifunction cryolipolysis machine. It has a central control panel with a touchscreen display. Several treatment heads are mounted on a rotating arm, each with a green light. The machine is set against a background of a light blue, textured surface.


360°
360° Refrigeration


Different sizes for changeable


Exchangeable cryo handles


Green light flashing


Adjustable time for cold & hot


Multifunction Cryolipolysis Machine

**360 ANGLE CRYOLIPOLYSIS
SLIMMING MACHINE**

Cryo handle screen



Home screen interface



360° ANGLE SURROUNDING CRYOLIPOLYSIS SLIMMING MACHINE

CRYO+Lipo Laser+ 40K cavitation
Body RF+Face RF



360° Refrigeration



Different sizes for
changeable



Adjustable time for
cold & hot



Exchangeable cryo
handles



Slimming



Green light flashing



- ❖ Cellulite removal
- ❖ Localized fat removal
- ❖ Lymph drained
- ❖ Skin tightening
- ❖ Body slimming, Reshape body line
- ❖ Pain relief for relaxation
- ❖ Improve blood circulation
- ❖ Combine cryolipolysis, cavitation treatment with RF to enhance slimming effect of beauty equipment

***** SMALL HANDLE DETAILS *****

Contours are securely fastened to the handle quickly,
very convenient for different areas treatment.




125*55*45mm



140*55*45mm



150*55*45mm



interchangeable

change contours with ease



Clasp

360° refrigeration



Treat area



arm



thigh



crus



CRIOlipOLYSIS



VACUUM



CAVITATION




LIPO LASER



RF



LED



4 HANDLES

Criolipolysis Handle

- * Minimum temperature of -16
- * Maximizes treatment area with 38%
- * Bottom switch control



Fast & Efficient



Freezing & Painless



Silent System



Size: 67*40mm
Unique double chin
treatment handle



Lipolaser handle



Radio Frequency
handles



Cavitation handle



40K Cavitation

direct into the fatty layer, speedily vibrate deep-seated cellulite, produce numberless vacuum cavitation, mightily strike the fatty cells, let them produce inner cracking, and dissolve to be the free fatty acid.



Body RF

uses proprietary technology to create a highly efficient and tightly woven energy matrix. This dense energy matrix penetrates multiple layers of the dermis heating it from the inside out. Complete skin penetration ensures exponentially improved clinical efficacy and products clearly visible results.



Face RF

RF radio frequency is one kind of high frequency electromagnetic wave. It can penetrate the dermis and subcutaneous tissue. Face RF can be used for face lifting, wrinkle removal.





Lipo laser

utilizes the latest low level/cold laser technology to reduce girth and spot fat over almost all areas of the body without surgery, downtime or redness. A typical 40-minute waistline treatment can reduce an individual by $\frac{1}{2}$ to $\frac{3}{4}$ of an inch while a full protocol of eight treatments usually results in multiple inch loss.



Fat on the abdomen



Thick arms




Thick legs



Fat on the waist


People who have excess fat on their abdomen, legs, arms, etc., who want to reshape beautiful lines without affecting their work and life.

How it works




STEP ONE

The 360° Cryo cup will target and cool fat cells to temperatures that trigger fat cell apoptosis.



STEP TWO

No damage to nerves or other tissues because lipids in fat crystallize at a warmer temperature than water in other cell types.



STEP THREE

allowing treatment, fat cells enter an apoptotic death sequence and are gradually removed in the next few weeks and months by the immune system.



STEP FOUR

Fat layer thickness significantly reduced.

real pictures




before and after



Customer feedback



Before



After



Before



After



Before



After



Before




After



Before



After



Before



After



Before



After



Before



After



Before



After

Cryo Technology

Cryolipolysis is a new, non-invasive way to gently and effectively reduce fat in targeted areas of the body that results in a noticeable, advanced-looking fat reduction in the treated areas. As triglyceride in fats will be converted into solid in particular low temperatures, it uses cooling technology to selectively target fat bulges and eliminate fat cells through a gradual process that does not harm the surrounding tissues, reduce unwanted fat. The contact cooling of the hand piece surface regulates the temperature of the skin and protects fine dermal structures, realizing the fast body-reshape effects while tightening skin!

Application:

- Body shaping: Shaping of the four limbs, buttocks, abdomen, recovery of Postpartum and liposuction surgery, Cellulite removal, Skin tightening.
- Body slimming: reduce the cellulite, reduce the fat, improve skin laxity.
- Lymph drained, Pain relief for relaxation, Improve blood circulation.
- Combine cryolipolysis, cavitation treatment with RF to enhance slimming effect of beauty equipment.

