


Alice Super Bubble Hydro Facial Machine

- 🌊 Bubble Cleaning
- 🌊 Skin Tighten
- 🌊 Fine Wrinkles
- 🌊 Deep Hydration
- 🌊 Moisturizing Nutrients
- 🌊 Improve Skin Texture



Nano Diamond Micro-carved Handle

Nano Diamond Micro-carved Handle



Scrubs away dead skin, softens cuticles, and promotes products absorption.

Painless and Non-invasive



Does not touch the dermis and will not damage the skin

Alice Super Bubble Machine

Deep Cleansing Skin System



Deep
Cleaning



Moisturizing
Nutrients Supply



Increase Skin
Transparency



Improve
Skin Texture



Firming
Skin



Diamond
Dermabrasion



Usage procedure diagram



1

Makeup Removal and Cleansing(2 Mins):

About 2 mins and then apply toner.



2

Water Bubble Cleaning (5-10 Mins):

Adjust the suitable strength of negative pressure and water output.



3

Diamond Dermabrasion (10-15 Mins):

Select the working head, mode and intensity of vibration.



4

High Pressure Spray (3-5 Mins):

Hyperbaric oxygen injection for enhanced water lock.



5

Carbon Oxygen Bubble (10-15 Mins):

Select the working head, mode and intensity of vibration.



6

Bipolar HR (8-10 Mins):

Regulating energy intensity and patterns.



7

Tripolar HR (8-10 Mins):

Regulating energy intensity and patterns.



8

Vibration Atomization (10-15 Mins):

Regulating energy intensity and pattern.



9

Hot and Cold (10-15 Mins):


Firming and contouring, redness and sensitivity relief.

Diamond Micro Dermabrasion Handle



Nanocrystals use nanoscale contacts that can be as thin as 250 nanometers and as short as 0.1 millimeters. Just open the stratum corneum of the skin without touching the dermis layer. It not only solves the problem of skin absorption, but also makes it painless and non-invasive!

Bubble Cleaning Handle



Vacuum negative pressure effectively removes impurities from different areas of the skin and moisturizes the skin multiple times.

Our Hydrofacial Machine

Non-invasive for Deeper Penetration and Longer Lasting Results



Steps One: CLEANSE

Exfoliate Skin Deep Clean for Effective Non-surgical Resurfacing Procedures.

Steps Two: EXOFOLIAT

Painless Gentle Suction to Remove the Particles Along with Any Dirt Debris, Imperfections, Blemishes, Wrinkles and Unwanted Pigmentation.

Steps Three: HYDRATE

Restore Cellular Activity to Achieve a Water-oil Balance, Regain Hydrated Skin to Promote Skin Health. Free Training and Online Video Support.

NewFace also offers products that work in conjunction with the corresponding serums, gels and masks to simultaneously brighten the skin and restore elasticity.

Specifications

Brand: NewFace

Voltage: 100~240V 50Hz 60Hz

Power: 100W

HR frequency: 0.2-0.5 Mhz

cold and hot scope: Cold 0-6C Hot 36-46C

Ultrasound: 1.1MHz

Display: 10.4 inches

Resolution: 1024X600

Clean water mill: Positive pressure 4-5Kg Negative pressure 96Kpa. Flow < 50L

Solution bottle



A. Deep layer cleaning

B. Whitening and moisturizing

C. Deep cleaning

D. Instrument piping cleaning

Probe: 3 soft heads, 8 hard heads




Diamond Hardhead : for pore thick skin

Silicone Soft Head: for sensitive mmskin

Click on the Suction: for oily skin

Welcome to inquire, online videos can be supported, free of charge.
Including how to use the different functions of the handle and the steps.

In the Box



Bubble Cleaning



Oxygen Bubble



Vibration
Atomization



Spray



Diamond
Dermabrasion



Cold & Hot



Tripolar HR



Bipolar HR

PRODUCT IMAGES

